210 Po

Nuclide Safety Data Sheet Polonium-210

radsafety@berkeley.edu

 $|^{210}$ Po

I. PHYSICAL DATA

Radiation¹: Gamma: 803 keV (<1%)

Alpha: 5,304 keV (100%)

Gamma Constant²:

Half-life [T_{1/2}]: Physical 138 days

Biological³: ~ 50 days; (excretion: 90% feces, 10% urine)

Effective: ~ 37 days

Specific Activity: 286,500 Ci/g [10,600 TBq/g] max.

Beta Range: Air: 610 cm [240 inches; 20 feet]

Water/Tissue: 0.76 cm [0.33 inches]
Plastic: 0.61 mm [3/8 inches]

II. RADIOLOGIAL DATA

Radiotoxicity⁴: 6.36E-12 Sv/Bq [0.023 mrem/uCi] of ¹⁴CO₂ inhaled;

5.64E-10 Sv/Bq [2.09 mrem/uCi] organic compounds inhaled/ingested Fat tissue [most labeled compounds]; bone [some labeled carbonates] Ingestion, inhalation, puncture, wound, skin contamination absorption

Radiological Hazard: External Exposure - None from weak ¹⁴C beta

Internal Exposure & Contamination - Primary concern

III. SHIELDING

Exposure Routes:

Critical Organ:

Shield ³²P with 3/8 inch Plexiglas and monitor for Bremstrahlung; If Bremstrahlung X-rays detected outside Plexiglas, apply 1/8 to 1/4 inch lead [Pb] shielding outside Plexiglas

IV. DOSIMETRY MONITORING

Always wear radiation dosimetry monitoring badges [body and/or ring] if required by your RUA.

V. DECTION & MEASUREMENT

Portable Survey Meters: Geiger-Mueller [e.g. Ludlum 44-9]

Wipe Test: Gamma Counter, Well Gamma Counter, or Liquid Scintillation Counter.

VI. SPECIAL PRECAUTIONS

- * Because of low gamma constant, no need to store ²¹⁰Po behind lead (Pb) shielding unless accessible dose rates >2 mR/hr observed.
- * Use tools to indirectly handle unshielded sources and potentially contaminated vessels; avoid direct hand contact.
- * Ensure that an appropriate, operational survey meter (e.g. Ludlum 44-9) is present in the work area and turned on whenever ²¹⁰Po is handled, so that any external exposure issues will be immediately apparent and hence quickly addressed.
- * Shield waste containers as needed to maintain accessible dose rate ALARA and < 2 mR/hr.</p>

VII. LAB PRACTICES

- 1. Disposable gloves, lab coats, and safety glasses are the minimum PPE [Personal Protective Equipment] required when handling radioactive material. Remove & discard potentially contaminated PPE prior to leaving the area where radioactive material is used.
- 2. Clearly outline radioactive material use areas with tape bearing the legend "radioactive." Cover lab bench tops where radioactive material will be handled with plastic-backed absorbent paper; change this covering periodically and whenever it's contaminated. Alternately cover benches with thick plastic sheeting (i.e. painter's drop cloth), periodically wipe it clean and replace it if torn.
- 3. Label each unattended radioactive material container with the radioactive symbol, isotope, activity, and, except for waste, the ICN [inventory control number]. Place containers too small for such labels in larger labeled containers.
- 4. Handle radioactive solutions in trays large enough to contain the material in the event of a spill.
- 5. Never eat, drink, smoke, handle contact lenses, apply cosmetics, or take/apply medicine in the lab; keep food, drinks, cosmetics, etc. out of the lab entirely. Do not pipette by mouth.

¹ Delacroix et al, <u>Radiation Protection Dosimetry - Radionuclide and Radiation Protection Data Handbook 2002</u> (Kent, England: Nuclear Technology Publishing, 2002), p. 146

² Shleien et al, Eds. <u>Handbook of Health Phys. & Rad. Health, 3rd ed</u>. (Baltimore, MD: Williams & Wilkins, 1998), p. 6-13

³ NRC: NUREG/CR-4884 <u>Interpretation of Bioassay Measurements</u> (Upton, New York: Brookhaven National Laboratories, 1987), p. B-788149

⁴ Federal Guidance Report No. 11 (Oak Ridge, TN; Oak Ridge National Laboratory, 1988), p. 175, 149

- 6.
- 7.
- Never store [human] food and beverage in refrigerators/freezers used for storing radioisotopes.

 Prevent skin contact with skin-absorbable solvents containing radioactive material.

 Fume hoods and biological safety cabinets for use with non-airborne radioactive material must be labeled "Caution" 8. Radioactive Material".